

TEST CHANGE

NOTIFICATION DATE: JUNE 2, 2014

EFFECTIVE DATE: JUNE 30, 2014

STRIATIONAL (STRIATED MUSCLE) ANTIBODIES, SERUM
TEST ID: STR

EXPLANATION: The Reference range and specimen stability are changing.

CURRENT REFERENCE VALUE	NEW REFERENCE VALUE
<1:60	<1:120

CURRENT SPECIMEN STABILITY			NEW SPECIMEN STABILITY		
Specimen	Temperature	Time	Specimen	Temperature	Time
Serum	Refrigerate	14 days	Serum	Refrigerate	28 days
	Ambient	72 hours		Ambient	72 hours
	Frozen			Frozen	28 days

QUESTIONS: Contact your Mayo Medical Laboratories' Regional Manager or
Steven Monson, MML Laboratory Technologist Resource Coordinator
Telephone: 800-533-1710