

TEST CHANGE

NOTIFICATION DATE: JUNE 2, 2014 **EFFECTIVE DATE:** JUNE 30, 2014

GLUTAMIC ACID DECARBOXYLASE (GAD65) ANTIBODY ASSAY, SERUM TEST ID: GD65S

EXPLANATION: Specimen stability updates.

CURRENT SPECIMEN STABILITY			NEW SPECIMEN STABILITY		
Specimen	Temperature	Time	Specimen	Temperature	Time
Serum	Refrigerate	14 days	Serum	Refrigerate	28 days
	Ambient	72 hours		Ambient	72 hours
	Frozen			Frozen	28 days

QUESTIONS: Contact your Mayo Medical Laboratories' Regional Manager or Steven Monson, MML Laboratory Technologist Resource Coordinator Telephone: 800-533-1710