



**MAYO**  
**Mayo Medical Laboratories**  
**1-800-533-1710**

**NEW TEST ANNOUNCEMENT**  
**REFERRAL**

**NOTIFICATION DATE:** June 30, 2014  
**EFFECTIVE DATE:** July 7, 2014

**ENDOMYSIAL (EMA) IGG ANTIBODY TITER**  
**Test ID: FEGAT**

**EXPLANATION:** The following new orderable, referred to IMMCO Diagnostics, Inc., will be available July 7, 2014. This test will replace Test ID FEMA

**REFERRAL LAB CODE:** 110

**PROFILE INFORMATION:**

Test ID	Reporting Name	Available Separately	Always Performed
FEMAS	Endomysial (EMA) Ab (IgG) screen	No	Yes
FEMAT	Endomysial (EMA) Ab (IgG) titer	No	No

**TEST ALGORITHM:** Endomysial (EMA) antibody (IgG) titer is performed at an additional charge if the screen is positive.

**METHODOLOGY:** Immunofluorescence (IFA)

**REFERENCE VALUES:** Negative: <1:2.5 titer  
 Positive: > or = 1:2.5 titer

**SPECIMEN REQUIREMENTS:** Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 3 mL of serum frozen in a plastic vial.

**SPECIMEN STABILITY INFORMATION:**

Specimen Type	Temperature	Time
Serum	Frozen (preferred)	365 days
	Ambient	30 days
	Refrigerated	30 days

**FEE:** \$ 87.00 an additional fee of \$31.00 if titer is performed

**CPT CODE:** 86255; 86256 if titer is performed

**DAY(S) SET UP:** Monday through Friday      **ANALYTIC TIME:** 2 days

**QUESTIONS:** Contact Mary Erath, MML Referrals Supervisor  
 Telephone: 800-533-1710