

Reporting Title: Endomysial (EMA) antibody (IgG)

Performing Location: Immco Diagnostics

Specimen Requirements:

Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 3 mL of serum frozen in a plastic vial.

| Specimen Type | Temperature | Time |
|---------------|--------------------|----------|
| Serum | Frozen (preferred) | 365 days |
| | Ambient | 30 days |
| | Refrigerated | 30 days |

Result Codes:

| Result ID | Reporting Name | Type | Unit | LOINC® |
|-----------|----------------------------------|--------------|-------|--------|
| Z4104 | Endomysial (EMA) Ab (IgG) screen | Alphanumeric | titer | |
| Z4105 | Endomysial (EMA) Ab (IgG) titer | Alphanumeric | titer | |

Components:

| Test ID | Reporting Name | CPT Units | CPT Code | Always Performed | Orderable Separately |
|---------|----------------------------------|-----------|----------|------------------|----------------------|
| FEMAS | Endomysial (EMA) Ab (IgG) screen | | | Yes | No |
| FEMAT | Endomysial (EMA) Ab (IgG) titer | | | No | No |

CPT Code Information:

86255

Reference Values:

Negative: <1:2.5 titer

Positive: > or = 1:2.5 titer

Test Performed by: IMMCO Diagnostics, Inc.
640 Ellicott Street
3rd Floor
Buffalo, NY 14203