

**Reporting Title:** Citalopram, S  
**Performing Location:** Rochester

**Specimen Requirements:**

Container/Tube: Red top

Specimen Volume: 0.5 mL

Collection Instructions: Spin down within 2 hours of draw.

Specimen Type	Temperature	Time
Serum Red	Refrigerated (preferred)	14 days
	Ambient	14 days
	Frozen	14 days

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
83730	Citalopram, S	Numeric	ng/mL	34635-3

**CPT Code:** 1 x 80299**Reference Values:**

50-110 ng/mL