

Reporting Title: Oxymorphone Confirmation, U
Performing Location: Rochester

Specimen Requirements:

Collection Container/Tube: Plastic urine container
Submission Container/Tube: Plastic, 60-mL urine bottle
Specimen Volume: 20 mL

Collection Instructions:

1. Collect a random urine specimen.
2. Submit 20 mL in 1 plastic bottle.
3. No preservative.

Additional Information:

1. No specimen substitutions.
2. No STATS are accepted for this procedure.
3. For situations where chain of custody is required, a Chain-of-Custody Kit (Supply T282) is available. For chain-of-custody information, see OXYCX / Oxycodone with Metabolite Confirmation, Chain of Custody, Urine
4. Additional drug panels and specific requests are available. Call Mayo Medical Laboratories at 800-533-1710 or 507-266-5700.
5. If urine creatinine is required or adulteration of the sample is suspected, the following test should be requested, ADULT / Adulterants Survey, Urine.
6. For additional information, please refer to ADULT / Adulterants Survey, Urine.

Specimen Type	Temperature	Time
Urine	Refrigerated (preferred)	14 days
	Frozen	14 days
	Ambient	72 hours

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
62622	Oxymorphone-by LC-MS/MS	Alphanumeric	ng/mL	In Process
36024	Interpretation	Alphanumeric		In Process

CPT Code Information:

83925



Reference Values:

Negative

Cutoff concentrations:

Oxymorphone: 100 ng/mL