

Test Definition: CDAU5

Drug Abuse Survey with Confirmation, Panel 5, Urine

Reporting Title: Confirmed Drug Abuse Panel 5, U

Performing Location: Rochester

Specimen Requirements:

Collection Container/Tube: Plastic urine container Submission Container/Tube: Plastic, 60-mL urine bottle

Specimen Volume: 30 mL Collection Instructions:

- 1. Collect a random urine specimen.
- 2. Submit 30 mL in 1 plastic bottle.
- 3. No preservative. Additional Information:
- No specimen substitutions.
- 2. No STATS are accepted for this procedure.
- 3. For situations where chain of custody is required, a Chain-of-Custody Kit (Supply T282) is available. For chain-of-custody information, see CDA5X / Drug Abuse Survey with Confirmation, Panel 5, Chain of Custody, Urine.
- 4. Additional drug panels and specific requests are available. Call Mayo Medical Laboratories at 800-533-1710 or 507-266-5700.
- 5. If urine creatinine is required or adulteration of the sample is suspected, the following test should be requested, ADULT / Adulterants Survey, Urine. For additional information, please refer to ADULT / Adulterants Survey, Urine.

| Specimen Type | Temperature | Time | |
|---------------|--------------------------|----------|--|
| Urine | Refrigerated (preferred) | 14 days | |
| | Frozen | 14 days | |
| | Ambient | 72 hours | |

Result Codes:

| Result ID | Reporting Name | Туре | Unit | LOINC® |
|-----------|---|---------|-------|---------|
| 2573 | Amphetamines | Numeric | ng/mL | 20410-7 |
| | Also used by tests: CDAU2, CDAU7, CDAU3, CDAU, CDAU1 | | | |
| 21652 | Cocaine | Numeric | ng/mL | 20519-5 |
| | Also used by tests: CDAU2, CDAU7, CDAU1, CDAU, CDAU3 | | | |
| 2577 | Opiates | Numeric | ng/mL | 19295-5 |
| | Also used by tests: CDAU, CDAU7, CDAU2, CDAU1, IDOAU, CDAU3 | | | |

Test Definition: CDAU5

Drug Abuse Survey with Confirmation, Panel 5, Urine

| Result ID | Reporting Name | Туре | Unit | LOINC® |
|-----------|---|---------|-------|---------|
| 2578 | Phencyclidine Also used by tests: CDAU2, CDAU3, CDAU7, IDOAU, CDAU, CDAU1 | Numeric | ng/mL | 19659-2 |
| 2664 | Tetrahydrocannabinol Also used by tests: CDAU3, CDAU7, CDAU | Numeric | ng/mL | 3530-3 |
| 20672 | Chain of Custody Also used by tests: IDOAU, CDAU1, CDAU7, CDAU3, CDAU2, CDAU | Numeric | | N/A |

CPT Code: 1 × 80104

Reflex Tests:

| Test ID | Reporting Name | CPT Units | CPT Code | Always Performed | Orderable Separately |
|---------|--------------------------------|-----------|----------|---------------------|-------------------------|
| AMPHU | U Amphetamines Confirmation, U | | 82145 | No | Yes |
| COKEU | Cocaine and metabolite Conf, U | 1 | 82520 | No | Yes |
| OPATU | Opiate Confirmation, U | | 83925 | No | Yes |
| PCPU | Phencyclidine Confirmation, U | 1 | 83992 | No | Yes |
| THCU | Carboxy-THC Confirmation, U | 1 | 80299 | No | Yes |

Result Codes for Reflex Tests:

| Test ID | Result ID | Reporting Name | Туре | Unit | LOINC® |
|---------|-----------|------------------------------------|--------------|-------|---------|
| AMPHU | 2934 | Amphetamine-by GC/MS | Alphanumeric | ng/mL | 19346-6 |
| AMPHU | 29278 | Phentermine-by GC/MS | Alphanumeric | ng/mL | 3958-6 |
| AMPHU | 2550 | Methamphetamine-by GC/MS | Alphanumeric | ng/mL | 3780-4 |
| AMPHU | 29279 | Pseudoephedrine/Ephedrine-by GC/MS | Alphanumeric | ng/mL | 26785-6 |
| AMPHU | 29280 | MDA (Ecstasy metabolite)-by GC/MS | Alphanumeric | ng/mL | 18355-8 |
| AMPHU | 29281 | MDMA (Ecstasy)-by GC/MS | Alphanumeric | ng/mL | 19570-1 |
| AMPHU | 21197 | Amphetamines Interpretation | Alphanumeric | | 69050-3 |
| AMPHU | 20641 | Chain of Custody | Numeric | | N/A |



Test Definition: CDAU5

Drug Abuse Survey with Confirmation, Panel 5, Urine

| Test ID | Result ID | Reporting Name | Туре | Unit | LOINC® |
|---------|-----------|------------------------------|--------------|-------|------------|
| COKEU | 2922 | Cocaine-by GC/MS | Alphanumeric | ng/mL | 3398-5 |
| COKEU | 2902 | Benzoylecgonine-by GC/MS | Alphanumeric | ng/mL | 3394-4 |
| COKEU | 20671 | Cocaine Interpretation | Alphanumeric | | 69050-3 |
| COKEU | 28080 | Chain of Custody | Numeric | | N/A |
| OPATU | 6543 | Codeine-by LC-MS/MS | Alphanumeric | ng/mL | 3508-9 |
| OPATU | 21177 | Hydrocodone-by LC-MS/MS | Alphanumeric | ng/mL | 3681-4 |
| OPATU | 21179 | Hydromorphone-by LC-MS/MS | Alphanumeric | ng/mL | 9835-0 |
| OPATU | 21178 | Oxycodone-by LC-MS/MS | Alphanumeric | ng/mL | 11246-6 |
| OPATU | 35023 | Oxymorphone-by LC-MS/MS | Alphanumeric | ng/mL | In Process |
| OPATU | 6542 | Morphine-by LC-MS/MS | Alphanumeric | ng/mL | 3831-5 |
| OPATU | 21180 | Opiates Interpretation | Alphanumeric | | 69050-3 |
| OPATU | 21181 | Chain of Custody | Numeric | | N/A |
| PCPU | 6673 | Phencyclidine-by GC/MS | Alphanumeric | ng/mL | 16254-5 |
| PCPU | 21192 | Phencyclidine Interpretation | Alphanumeric | | 69050-3 |
| PCPU | 21193 | Chain of Custody | Numeric | | N/A |
| THCU | 2497 | Carboxy-THC- by GC/MS | Alphanumeric | ng/mL | 3530-3 |
| THCU | 21186 | Carboxy-THC Interpretation | Alphanumeric | | 69050-3 |
| THCU | 21187 | Chain of Custody | Numeric | | N/A |

Reference Values:

Negative

Screening cutoff concentrations Amphetamines: 500 ng/mL

Cocaine (benzoylecgonine-cocaine metabolite): 150 ng/mL

Opiates: 300 ng/mL Phencyclidine: 25 ng/mL

Tetrahydrocannabinol carboxylic acid: 20 ng/mL

This report is intended for use in clinical monitoring or management of patients. It is not intended for use in

employment-related testing.