
Reporting Title: 6-MAM Confirmation, U
Performing Location: Rochester**Specimen Requirements:**

Container/Tube: Plastic, 60-mL urine bottle

Specimen Volume: 20 mL

Collection Instructions:

1. Collect a random urine specimen.
2. No preservative.

Additional Information:

1. If urine creatinine is required or adulteration of the sample is suspected, the following test should be requested, ADULT / Adulterants Survey, Urine. For additional information, please refer to ADULT / Adulterants Survey, Urine.
2. For situations where chain-of-custody is required, a Chain-of-Custody Kit (Supply T282) is available. For chain-of-custody information, see 6MAMX / 6-Monoacetylmorphine (6-MAM), Chain of Custody, Urine.

Specimen Type	Temperature	Time
Urine	Refrigerated (preferred)	14 days
	Frozen	14 days
	Ambient	72 hours

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
30983	6-Monoacetylmorphine by GC/MS	Alphanumeric	ng/mL	10975-1
29898	6-MAM Interpretation	Alphanumeric		59462-2
29899	Chain of Custody	Numeric		N/A

CPT Code: 1 x 83925**Reference Values:**

Negative

Cutoff concentrations:

6-MAM

<5 ng/mL