



**Reporting Title:** Confirmed Drug Abuse Panel9, CoC, U

**Performing Location:** Rochester

**Specimen Requirements:**

Container/Tube: Chain-of-Custody Kit (Supply T282) containing the specimen containers, seals, and documentation required

Specimen Volume: 30 mL

Collection Instructions: Collect specimen in the container provided, seal, and submit with the associated documentation to satisfy the legal requirements for chain-of-custody testing.

Specimen Type	Temperature	Time
Urine	Refrigerated (preferred)	14 days
	Frozen	14 days
	Ambient	72 hours

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
36262	Alcohol	Numeric	mg/dL	In Process
36253	Amphetamines	Numeric	ng/mL	In Process
36258	Barbiturates	Numeric	ng/mL	In Process
36259	Benzodiazepines	Numeric	ng/mL	In Process
36254	Cocaine	Numeric	ng/mL	In Process
36260	Methadone	Numeric	ng/mL	In Process
36255	Opiates	Numeric	ng/mL	In Process
36256	Phencyclidine	Numeric	ng/mL	In Process
36257	Tetrahydrocannabinol	Numeric	ng/mL	In Process
36261	Chain of Custody	Numeric		N/A

**CPT Code:** 1 x 80104



**Reflex Tests:**

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
OPATX	Opiate Confirmation, CoC, U	1	83925	No	Yes
COCH	Chain of Custody Processing	1	99001	Yes	No
ADLTX	Adulterants Survey, CoC, U	1	81005	Yes	Yes
AMPHX	Amphetamines Confirmation, CoC, U	1	82145	No	Yes
BARBX	Barbiturates Confirmation, CoC, U	1	82205	No	Yes
BENZX	Benzodiazepines Conf, CoC, U	1	80154	No	Yes
COKEX	Cocaine and metabolite Conf, CoC, U	1	82520	No	Yes
ETOHX	Ethanol, CoC, U	1	82055	No	Yes
MTDNX	Methadone Confirmation, CoC, U	1	83840	No	Yes
PCPX	Phencyclidine Confirmation, CoC, U	1	83992	No	Yes
THCX	Carboxy-THC Confirmation, CoC, U	1	80299	No	Yes

**Result Codes for Reflex Tests:**

Test ID	Result ID	Reporting Name	Type	Unit	LOINC®
OPATX	6541	Opiates Immunoassay Screen	Numeric		19295-5
OPATX	36213	Codeine-by LC-MS/MS	Alphanumeric	ng/mL	In Process
OPATX	36214	Hydrocodone-by LC-MS/MS	Alphanumeric	ng/mL	In Process
OPATX	36215	Hydromorphone-by LC-MS/MS	Alphanumeric	ng/mL	In Process
OPATX	36216	Oxycodone-by LC-MS/MS	Alphanumeric	ng/mL	In Process
OPATX	36212	Oxymorphone-by LC-MS/MS	Alphanumeric	ng/mL	In Process
OPATX	36217	Morphine-by LC-MS/MS	Alphanumeric	ng/mL	In Process
OPATX	36218	Opiates Interpretation	Alphanumeric		In Process
OPATX	36219	Chain of Custody	Numeric		N/A
ADLTX	36121	Creatinine, U	Numeric	mg/dL	In Process
ADLTX	36122	Specific Gravity	Numeric		In Process
ADLTX	36123	pH	Numeric		In Process
ADLTX	36124	Oxidants	Alphanumeric		In Process



Test ID	Result ID	Reporting Name	Type	Unit	LOINC®
ADLTX	36125	Nitrites	Alphanumeric	mcg/mL	In Process
ADLTX	36126	Comment	Alphanumeric		In Process
AMPHX	6538	Amphetamines Immunoassay Screen	Numeric		19343-3
AMPHX	36128	Amphetamine-by GC/MS	Alphanumeric	ng/mL	In Process
AMPHX	36129	Phentermine-by GC/MS	Alphanumeric	ng/mL	In Process
AMPHX	36130	Methamphetamine-by GC/MS	Alphanumeric	ng/mL	In Process
AMPHX	36131	Pseudoephedrine/Ephedrine-by GC/MS	Alphanumeric	ng/mL	In Process
AMPHX	36132	MDA (Ecstasy metabolite)-by GC/MS	Alphanumeric	ng/mL	In Process
AMPHX	36133	MDMA (Ecstasy)-by GC/MS	Alphanumeric	ng/mL	In Process
AMPHX	36134	Amphetamines Interpretation	Alphanumeric		In Process
AMPHX	36135	Chain of Custody	Numeric		N/A
BARBX	6668	Barbiturates Immunoassay Screen	Numeric		19270-8
BARBX	36143	Butalbital-by GC/MS	Alphanumeric	ng/mL	In Process
BARBX	36144	Amobarbital-by GC/MS	Alphanumeric	ng/mL	In Process
BARBX	36145	Pentobarbital-by GC/MS	Alphanumeric	ng/mL	In Process
BARBX	36146	Secobarbital-by GC/MS	Alphanumeric	ng/mL	In Process
BARBX	36147	Phenobarbital-by GC/MS	Alphanumeric	ng/mL	In Process
BARBX	36148	Barbiturates Interpretation	Alphanumeric		In Process
BARBX	36149	Chain of Custody	Numeric		N/A
BENZX	6670	Benzodiazepines Immunoassay Screen	Numeric		14316-4
BENZX	36150	Nordiazepam-by GC/MS	Alphanumeric	ng/mL	In Process
BENZX	36151	Oxazepam-by GC/MS	Alphanumeric	ng/mL	In Process
BENZX	36152	Lorazepam-by GC/MS	Alphanumeric	ng/mL	In Process
BENZX	36153	Temazepam-by GC/MS	Alphanumeric	ng/mL	In Process
BENZX	36154	OH-Ethyl-Flurazepam-by GC/MS	Alphanumeric	ng/mL	In Process
BENZX	36155	7-NH-Clonazepam-by GC/MS	Alphanumeric	ng/mL	In Process
BENZX	36156	Alpha OH-Alprazolam-by GC/MS	Alphanumeric	ng/mL	In Process
BENZX	36157	7-NH-Flunitrazepam-by GC/MS	Alphanumeric	ng/mL	In Process
BENZX	36158	Alpha OH-Triazolam-by GC/MS	Alphanumeric	ng/mL	In Process
BENZX	36159	Benzodiazepines Interpretation	Alphanumeric		In Process
BENZX	36160	Chain of Custody	Numeric		N/A



Test ID	Result ID	Reporting Name	Type	Unit	LOINC®
COKEX	2903	Cocaine Immunoassay Screen	Numeric		3393-6
COKEX	36162	Cocaine-by GC/MS	Alphanumeric	ng/mL	In Process
COKEX	36163	Benzoyllecgonine-by GC/MS	Alphanumeric	ng/mL	In Process
COKEX	36164	Cocaine Interpretation	Alphanumeric		In Process
COKEX	36165	Chain of Custody	Numeric		N/A
ETOHX	30927	Ethanol Immunoassay Screen	Numeric	mg/dL	5645-7
ETOHX	36187	Ethanol-by GC-FID	Numeric	mg/dL	In Process
ETOHX	36188	Ethanol Interpretation	Alphanumeric		In Process
ETOHX	36189	Chain of Custody	Numeric		N/A
MTDNX	21105	Methadone Immunoassay Screen	Numeric		16199-2
MTDNX	36208	EDDP-by GC-MS	Alphanumeric	ng/mL	In Process
MTDNX	36209	Methadone-by GC-MS	Alphanumeric	ng/mL	In Process
MTDNX	36210	Methadone Interpretation	Alphanumeric		In Process
MTDNX	36211	Chain of Custody	Numeric		N/A
PCPX	6672	Phencyclidine Immunoassay Screen	Numeric		19659-2
PCPX	36229	Phencyclidine-by GC/MS	Alphanumeric	ng/mL	In Process
PCPX	36230	Phencyclidine Interpretation	Alphanumeric		In Process
PCPX	36231	Chain of Custody	Numeric		N/A
THCX	2449	Carboxy-THC Immunoassay Screen	Numeric		19415-9
THCX	36238	Carboxy-THC- by GC/MS	Alphanumeric	ng/mL	In Process
THCX	36239	Carboxy-THC Interpretation	Alphanumeric		In Process
THCX	36240	Chain of Custody	Numeric		N/A

**Reference Values:**

Negative

Screening cutoff concentrations

Amphetamines: 500 ng/mL

Barbiturates: 200 ng/mL

Benzodiazepines: 200 ng/mL

Cocaine (benzoyllecgonine-cocaine metabolite): 150 ng/mL

Ethanol: 30 mg/dL

Methadone: 150 ng/mL

Opiates: 300 ng/mL

Phencyclidine: 25 ng/mL

Tetrahydrocannabinol carboxylic acid: 20 ng/mL



---

This report is intended for use in clinical monitoring or management of patients. It is not intended for use in employment-related testing.