

Patient Name REPORTVALIDATION,AUTOMATION D...	Patient ID RVDMOHB036	Age 40	Gender F	Order # RVDMOHB036
Ordering Phys				DOB 01/01/1971
Client Order # RVDMOHB036	Account Information C7028846-DLMP Rochester SDSC 2 - Client Support Rochester, MN 55901			Report Notes
Collected 11/14/2011 16:22				
Printed 10/21/2013 13:02				

Test	Flag	Results	Unit	Reference Value	Perform Site*
------	------	---------	------	-----------------	---------------

BCR/ABL, p190, Quant, Monitor
RECEIVED: 11/15/2011 07:56 **REPORTED:** 11/15/2011 07:56

Specimen Type

Peripheral blood

MCR

Final Diagnosis:

MCR

Peripheral blood, BCR/ABL mRNA level analysis (p190 fusion form):

Positive. BCR/ABL p190 mRNA transcripts were detected and estimated to represent 1.20% of total abl (%bcr/abl(p190):abl).

Previous specimens from this patient have been tested in this laboratory using this method and a summary of results will be faxed provided our laboratory has a current, secure fax number on file. If a faxed report is not received, please call the Mayo Clinic Molecular Hematopathology Laboratory (1-507-266-0489) to provide a fax number.

Signing Pathologist: Carey Lueck

Method summary-BCR/ABL, p190 fusion: BCR/ABL p190 mRNA transcript level was evaluated using quantitative, reverse transcription PCR. The assay detects the most common fusion form (e1/a2), but does not detect other fusions, including the p210, which is the most common form found in chronic myelogenous leukemia. This assay should only be ordered for monitoring patients with a previously identified p190 fusion form. Test #89006(BCR/ABL, mRNA detection, RT-PCR, Quantitative, Diagnostic) should be ordered if the test is being performed in a diagnostic setting and test #89007(BCR/ABL, p210, mRNA detection, RT-PCR, Quantitative, Monitoring CML) should be ordered if this patient is being monitored for a known p210 fusion form. See the Mayo Medical Laboratories Interpretive Handbook for method details. Please contact the Mayo Molecular Hematopathology Laboratory at 507-266-0489 with questions or if additional testing is required.

Typical clinical samples have detection limits ranging from 0.01% to 0.0001% bcr/abl:abl. Most patients at diagnosis have a bcr/abl:abl result in the range of 20% to 100%.

The reproducibility of this assay is such that results within 0.5 log should be considered equivalent. Trends in the level of BCR/ABL mRNA should be followed and clinically

Performing Site Legend on Last Page of Report

Patient Name REPORTVALIDATION,AUTOMATION D...	Collection Date and Time 11/14/2011 16:22	Report Status Final
Page 1 of 2		>> Continued on Next Page >>

* Report times for Mayo performed tests are CST/CDT

Patient Name REPORTVALIDATION,AUTOMATION D...	Patient ID RVDMOHB036	Age 40	Gender F	Order # RVDMOHB036
Ordering Phys				DOB 01/01/1971
Client Order # RVDMOHB036	Account Information			Report Notes
Collected 11/14/2011 16:22	C7028846-DLMP Rochester SDSC 2 - Client Support			
Printed 10/21/2013 13:02	Rochester, MN 55901			

Test	Flag	Results	Unit	Reference Value	Perform Site*
------	------	---------	------	-----------------	---------------

significant changes verified with a subsequent specimen.
 Laboratory developed test.

* Performing Site:

MCR	Mayo Clinic Dpt of Lab Med & Pathology 200 First St SW Rochester, MN 55905	Lab Director:
-----	---	---------------

Patient Name REPORTVALIDATION,AUTOMATION D...	Collection Date and Time 11/14/2011 16:22	Report Status Final
Page 2 of 2		** End of Report **

* Report times for Mayo performed tests are CST/CDT