

**ANTI-IGA**  
**#91993**

**EXPLANATION OF CHANGE:** Effective July 2, 2010 Unit Code 91993, referred to IBT Reference Laboratory, will become obsolete.

**ALTERNATIVE MAYO TEST:** #8154, Anti-IgA

**METHODOLOGY:** Fluorescent Microsphere Immunoassay

**SPECIMEN REQUIREMENTS:** Draw blood in a plain, red-top tube(s) or a serum gel tube(s). Spin down and send 1 mL of serum frozen in plastic vial.

**Note:** If testing is due to a transfusion reaction, the specimen should be drawn approximately a minimum of 10 days after the reaction. Testing performed on serum drawn immediately after a reaction may be falsely negative. Transfused IgA may deplete anti-IgA antibodies.

**LIST FEE:** \$138.90

**CPT CODE:** 83520

**DAY(S) SET UP:** Tuesday; Thursday

**ANALYTIC TIME:** 1 day

**QUESTIONS:** Contact Mayo Medical Laboratories  
Telephone: 800-533-1710