

**C1 ESTERASE INHIBITOR ANTIGEN, SERUM
#200132**

NOTE: The testing and reporting is equivalent to Rochester test 8198, C1 Esterase Inhibitor Antigen, S.

USEFUL FOR:

- Diagnosis of hereditary angioedema
- Monitoring levels of C1 esterase inhibitor in response to therapy

METHODOLOGY: Nephelometry

REFERENCE VALUES: 19-37 mg/dL

SPECIMEN REQUIREMENTS: Draw blood in a plain, red-top tube(s). Fasting preferred but not required. Spin down and send 1 mL of serum frozen in plastic vial.

CAUTIONS: Quantitation of specific proteins by nephelometric means may not be possible in lipemic sera due to the extreme light scattering properties of the specimen. Turbidity and particles in the specimen may result in extraneous light scattering signals, resulting in variable specimen analysis.

LIST FEE: \$128.60

CPT CODE: 83520

ANALYTIC TIME: 1 Day

DAY(S) SET-UP: Mon-Fri, 3rd Shift

QUESTIONS: Contact your Mayo Medical Laboratories' Regional Manager or
Mary Erath, MML Laboratory Technologist Resource Coordinator
Telephone: 800-533-1710