

PIVKA-II
#91493

EXPLANATION OF CHANGE: Effective immediately, Unit code 91493, referred to Esoterix Coagulation, will be updated to reflect the following. The normal range change is due to variation of the reagent lot to lot performance.

	CURRENT	NEW
Reference Range	18y: Use published ranges 18-60y: <3.6 **RUO-When making a clinical diagnosis, this test should be used with other established laboratory tests or clinical observations.	<18y: Use published ranges 18-60y: <6.3 ng/mL **RUO-When making a clinical diagnosis, this test should be used with other established laboratory tests or clinical observations.

QUESTIONS: Contact Mayo Medical Laboratories
 Telephone: 800-533-1710