



MAYO CLINIC
 Mayo Medical Laboratories

TEST DEFINITION

03/11/2011

MML Name/order code cross-reference

Code Name

 89379 Thyroglobulin, FNA, Lymph Node

MML MML Test setup information

ORDER CODE	EFF DATE	TC	TITLE	Checking Normals	Print normals (# coded)	Perform Site *
89379	3/15/2011		Thyroglobulin, FNA, Lymph Node			SDL
			Transport temp : Frozen\Refrig OK\Ambient NO			
			89379 Thyroglobulin, FNA, Lymph Node			
			Units: ng/mL			
			NO SEX			
			All Ages :		; <=1.0	
			MALE			
			All Ages :		; <=1.0	
			FEMALE			
			All Ages :		; <=1.0	

TEST CODE ALWAYS MESSAGE - [Z89379]
 Z89379 Tumor marker tests are not specific for malignancy. Correlation of this test result with imaging and cytology is recommended. This test was developed and its performance characteristics determined by Laboratory Medicine and Pathology, Mayo Clinic, Rochester, MN. This test has not been cleared or approved by the U.S. Food and Drug Administration. The testing method is an immunoenzymatic assay manufactured by Beckman Coulter Inc. and performed on the UniCel DxI 800. Values obtained with different assay methods or kits may be different and cannot be used interchangeably.

MML MML Test setup information

*Performing Site Legend

=====
 SDL Mayo Clinic Dpt of Lab Med & Pathology LAB DIRECTOR: Franklin R. Cockerill, III, M.D.
 Superior Dr.
 3050 Superior Dr. NW
 Rochester, MN 55901
 =====

MML
CODE

TEXT

Messages used as normals

Total of 0 normals codes

*** End of Report ***