



MAYO
Mayo Medical Laboratories
1-800-533-1710

NEW TEST ANNOUNCEMENT
NOTIFICATION DATE: November 22, 2011
EFFECTIVE DATE: December 15, 2011

EVEROLIMUS, BLOOD
Test ID: EVERO
Secondary ID: 60537

USEFUL FOR: Management of everolimus immunosuppression in solid organ transplant.

METHODOLOGY: Liquid Chromatography-Tandem Mass Spectrometry (LC-MS/MS)

REFERENCE VALUES: 3-8 ng/mL

SPECIMEN REQUIREMENTS:

Container/Tube: Lavender top (EDTA)

Specimen Volume: 3 mL

Collection Instructions:

1. Draw blood immediately before a schedule dose.
2. **Do not centrifuge.**
3. Send specimen in original tube.

Additional Information: Therapeutic range applies to trough specimens drawn immediately prior to a.m. dose.

SPECIMEN STABILITY:

Specimen Type	Temperature	Time
Whole Blood EDTA	Refrigerated (preferred)	14 days
	Ambient	12 hours
	Frozen	14 days

LIST FEE: \$206.80

CPT CODE: 80299

DAY(S) SET UP: Monday through Sunday 1 p.m. **ANALYTIC TIME:** Same Day/1 day

QUESTIONS: Contact your Mayo Medical Laboratories' Regional Manager or
Jim Nielsen, MML Laboratory Technologist Resource Coordinator
Telephone: 800-533-1710