

## SPECIMEN STABILITY/FILE DEFINITION CHANGE

MML JACKSONVILLE

NOTIFICATION DATE: October 23, 2012 EFFECTIVE DATE: November 16, 2012

## SIROLIMUS, BLOOD

Test ID: 800072

**EXPLANATION:** There will be a change in the specimen stability information for Sirolimus, Blood. Refrigerate specimens will be changing from 14 days to 7 days.

## **NEW SPECIMEN STABILITY INFORMATION:**

Specimen Type	Temperature	Time
Whole Blood EDTA	Refrigerated (preferred)	7 days
	Ambient	12 hours

QUESTIONS: Contact your Mayo Medical Laboratories' Regional Manager or Greg Renkly, MML Laboratory Technologist Resource Coordinator Telephone: 800-533-1710