



MAYO
Mayo Medical Laboratories
1-800-533-1710

NEW TEST ANNOUNCEMENT
NOTIFICATION DATE: December 4, 2012
EFFECTIVE DATE: January 15, 2013

IOHEXOL, PLASMA

Test ID: HEXP

USEFUL FOR: Determining glomerular filtration (GFR)

NOTE: Testing is for approved accounts only

METHODOLOGY: Liquid Chromatography-Tandem Mass Spectrometry (LC-MS/MS)

REFERENCE VALUES: Not applicable

SPECIMEN REQUIREMENTS:

Container/Tube: Plastic, 6-mL tube (Supply T465)

Specimen Volume: 1 mL plasma from timed specimens collected in Sodium Heparin tube

SPECIMEN STABILITY INFORMATION

Specimen Type	Temperature	Time
Plasma Heparin	Refrigerated (preferred)	7 days
	Ambient	Not acceptable
	Frozen	35 days

LIST FEE: Please contact your Regional Manager for your account's fee information

CPT CODES:
83789-LC/MS/MS

DAY(S) SET UP: Monday through Friday; 8 a.m.-4:30 p.m. **ANALYTIC TIME:** 3 days

QUESTIONS: Contact your Mayo Medical Laboratories' Regional Manager or
Richard R. Einerson, MML Laboratory Technologist Resource Coordinator
Telephone: 800-533-1710