
Reporting Title: Fusarium oxysporum/vasinfectum IgE
Performing Location: VIRACOR-IBT LABORATORIES, INTERFACE

Specimen Requirements:

Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

Specimen Type	Temperature	Time
Serum	Refrigerated (preferred)	28 days
	Frozen	365 days
	Ambient	28 days

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
Z3090	Fusarium oxysporum/vasinfectum IgE	Alphanumeric		In Process
Z3091	Class	Alphanumeric		In Process

CPT Code Information:

86003

Reference Values:

<0.35 kU/L

Test Performed by: Viracor-IBT Laboratories
1001 NW Technology Dr.
Lee's Summit, MO 64086