
Reporting Title: JAK2 Exon 12 Mutation Detection, B
Performing Location: Rochester

Specimen Requirements:

Specimen must arrive within 72 hours of draw.

Container/Tube:

Preferred: Lavender top (EDTA)

Acceptable: ACD

Specimen Volume: 4 mL

Collection Instructions:

1. Invert several times to mix blood.
2. Send specimen in original tube.

Additional Information: Date of draw is required.

Forms: Hematopathology Patient Information Sheet (Supply T676) in Special Instructions

Specimen Type	Temperature	Time
Whole blood	Ambient (preferred)	72 hours
	Refrigerated	72 hours

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
20194	Final Diagnosis:	Alphanumeric		34574-4

CPT Code: 1 x 81403

Reference Values:

An interpretive report will be provided.