

Reporting Title: Pepper Bell/Paprika IgG

Performing Location: VIRACOR-IBT LABORATORIES, INTERFACE

Specimen Requirements:

Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

| Specimen Type | Temperature | Time |
|---------------|--------------------------|----------|
| Serum | Refrigerated (preferred) | 28 days |
| | Frozen | 365 days |
| | Ambient | 7 days |

Result Codes:

| Result ID | Reporting Name | Type | Unit | LOINC® |
|-----------|-------------------------|--------------|------|------------|
| FPBPG | Pepper Bell/Paprika IgG | Alphanumeric | | In Process |

CPT Code Information:

86001

Reference Values:

<2 mcg/mL

The reference range listed on the report is the lower limit of quantitation for the assay. The clinical utility of food-specific IgG tests has not been established. These tests can be used in special clinical situations to select foods for evaluation by diet elimination and challenge in patients who have food-related complaints. It should be recognized that the presence of food-specific IgG alone cannot be taken as evidence of food allergy and only indicates immunologic sensitization by the food allergen in question. This test should only be ordered by physicians who recognize the limitations of the test.

Test Performed by: Viracor-IBT Laboratories
1001 NW Technology Dr.
Lee's Summit, MO 64086