

Reporting Title: HLAB 1502 Carbamazepine, Saliva**Performing Location:** Rochester**Specimen Requirements:**

Multiple saliva genotype tests can be performed on a single specimen after a single extraction. See Multiple Saliva Genotype Tests in Special Instructions for a list of tests that can be ordered together.

Container/Tube: Oragene DNA Self-Collection Kit (Supply T651)

Specimen Volume: Full tube

Collection Instructions:

1. Fill tube to line.
2. Send specimen in original container per kit instructions.

Forms: New York Clients-Informed consent is required. Please document on the request form or electronic order that a copy is on file. An Informed Consent for Genetic Testing (Supply T576) is available in Special Instructions.

Specimen Type	Temperature	Time
Saliva	Ambient	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
33052	HLA-B 1502 Result	Alphanumeric		In Process
33053	HLA-B 1502 Reviewed by	Numeric		N/A
33054	HLA-B 1502 Interpretation	Alphanumeric		69047-9

CPT Code: 1 x 81381**Reference Values:**

An interpretive report will be provided.