
Reporting Title: C2 Complement, Functional, S, NR
Performing Location: Rochester

Specimen Requirements:

Collection Container/Tube: Red top

Submission Container/Tube: Plastic vial

Specimen Volume: 1 mL

Collection Instructions:

1. Immediately after drawing the specimen, place the tube on wet ice.
2. Spin down and separate serum from clot.
3. Immediately freeze specimen.

Additional Information: Fasting preferred.

Specimen Type	Temperature	Time
Serum Red	Frozen	21 days

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
C2FX	C2 Complement,Functional,S Also used by tests: C2	Numeric	U/mL	33408-6
INT53	Interpretation Also used by tests: C2	Alphanumeric		N/A

CPT Code: 1 x 86161**Reference Values:**

25-47 U/mL