
Reporting Title: HBs Antigen Prenatal, S
Performing Location: Jacksonville

Specimen Requirements:

Collection Container/Tube: Serum gel

Submission Container/Tube: Plastic vial

Specimen Volume: 2 mL

Collection Instructions: Spin down and remove serum from clot within 24 hours.

Additional Information: Date of draw is required.

Specimen Type	Temperature	Time
Serum	Frozen (preferred)	
	Refrigerated	7 days
	Ambient	24 hours

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
800345	HBs Antigen Prenatal, S	Alphanumeric		In Process

Reference Values:

Negative

See Viral Hepatitis Serologic Profiles in Special Instructions.