

**Reporting Title:** Euglobulin Clot Lysis Time

**Performing Location:** Quest SJC, Interface

**Specimen Requirements:**

Collect blood in 3.2% sodium citrate light blue top tube. Centrifuge within 30 minutes after collection, freeze immediately. Send 2 mL platelet-poor plasma in plastic vial frozen.

Note:

1. Prohibit exercise prior to drawing sample.
2. To avoid release of plasminogen activator, do not massage vein vigorously, pump fist excessively or leave tourniquet in place for a prolonged period.

Specimen Type	Temperature	Time
Plasma Na Cit	Frozen	21 days

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
FECLT	Euglobulin Clot Lysis Time	Alphanumeric		3244-1

**CPT Code Information:**

85360

**Reference Values:**

>60 min

Test Performed by: Quest Diagnostics Nichols Institute  
33608 Ortega Highway  
San Juan Capistrano, CA 92690