

**TEST OBSOLETE****NOTIFICATION DATE:** May 20, 2014**EFFECTIVE DATE:** Immediately**C-PEPTIDE, URINE**  
Test ID: FCPEP

**EXPLANATION:** Per notification received from Esoterix Endocrinology, effective immediately, Test ID FCPEP, will become obsolete.

**RECOMMENDED ALTERNATIVE TEST:**

ZW76-Miscellaneous LabCorp Burlington

Test name: C-Peptide, 24-Hour Urine

Referral lab code: 003236

**METHODOLOGY:** Electrochemiluminescence immunoassay (ECLIA)**SPECIMEN REQUIREMENTS:**

Instruct the patient to void at 8 AM and discard the specimen. Then collect all urine including the final specimen voided at the end of the 24-hour collection period (ie, 8 AM the next morning) in a 24-hour urine collection container. pH should be 4-7. Ship a 10 mL aliquot of the 24-hour urine collection, frozen in a plastic container.

Note: In patients receiving therapy with high biotin doses (ie, >5 mg/day), no sample should be taken until at least eight hours after the last biotin administration.

**SPECIMEN STABILITY INFORMATION:**

Specimen Type	Temperature	Time
Urine	Frozen (preferred)	14 days
	Refrigerated	14 days
	Ambient	1 day

**FEE:** \$50.00**CPT CODE:** 84681**DAY(S) SET UP:** Monday through Friday**ANALYTIC TIME:** 2 days