

**TEST STATUS**

NOTIFICATION DATE: July 16, 2014

EFFECTIVE DATE: August 25, 2014

**OXYCODONE CONFIRMATION, QUANTITATIVE, URINE**

Test ID: FOXYC

**EXPLANATION:** Effective August 25, 2014, Test ID FOXYC, referred to Medtox Laboratories, will become obsolete.

**RECOMMENDED ALTERNATIVE TEST:****OXYCODONE WITH METABOLITE CONFIRMATION, URINE**

Test ID: OXYCU

**USEFUL FOR:** Detection and quantification of oxycodone and oxymorphone in urine.

**METHODOLOGY:** Liquid Chromatography-Tandem Mass Spectrometry (LC-MS/MS)

**REFERENCE VALUES:** Negative

Cutoff concentrations:

Oxycodone: 100 ng/mL

Oxymorphone: 100 ng/mL

**SPECIMEN REQUIREMENTS:**

**Collection Container/Tube:** Plastic urine container

**Submission Container/Tube:** Plastic, 60-mL urine bottle (Supply T313)

**Specimen Volume:** 20 mL

**Minimum Volume:** 2.5 mL

**Collection Instructions:**

1. Collect a random urine specimen.
2. Submit 20 mL in 1 plastic bottle.
3. No preservative.

**Additional Information:**

1. No specimen substitutions.
2. No STATS are acceptable for this procedure.
3. For situations where chain of custody is required, a Chain-of-Custody Kit (Supply T282) is available. For chain-of-custody information, see OXYCX / Oxycodone with Metabolite Confirmation, Chain of Custody, Urine.
4. Additional drug panels and specific requests are available. Call Mayo Medical Laboratories at 800-533-1710.
5. If urine creatinine is required or adulteration of the sample is suspected, the following test should be requested, ADULT / Adulterant Survey, Urine.
6. For additional information, please refer to ADULT / Adulterant Survey, Urine.

**SPECIMEN STABILITY INFORMATION:**

Specimen Type	Temperature	Time
Urine	Refrigerated (preferred)	14 days
	Ambient	72 hours
	Frozen	14 days

**CAUTIONS:** Other drugs in the opioid class, such as fentanyl, meperidine, methadone, and opiate antagonists such as naloxone, are not detected.

**CPT CODE:** 83925

**DAYS SET UP:** Monday through Friday

**ANALYTIC TIME:** 2 days

**QUESTIONS:** Contact Mary Erath, MML Referrals Supervisor  
Telephone: 800-533-1710