

**SPECIMEN REQUIREMENT
UPDATE**

NOTIFICATION DATE: June 2, 2009

EFFECTIVE DATE: July 1, 2009

**Lysozyme (Muramidase), Urine
#8417**

EXPLANATION: Random urine collections will no longer be an acceptable alternative specimen type for #8417, Lysozyme (Muramidase), Urine.

Specimen Required:

12 mL from a 24-hour urine collection. No preservative. Refrigerate specimen during collection, and send specimen frozen in a plastic, 13-mL urine tube.

Note: 24-Hour volume is required on request form for processing.