

SPECIMEN REQUIREMENT CHANGE

NOTIFICATION DATE: June 16, 2010 **EFFECTIVE DATE:** Immediately

LACTATE DEHYDROGENASE (LD), SERUM #8344

LACTATE DEHYDROGENASE (LD) ISOENZYMES, SERUM #8679

EXPLANATION: The preferred transport temperature for #8344, Lactate Dehydrogenase (LD), will change. Effective immediately, the specimen transport temperature requirements for #8344, Lactate Dehydrogenase (LD), Serum, and #8679, Lactate Dehydrogenase (LD) Isoenzymes, Serum, will be the same.

CURRENT PREFERRED TRANSPORT TEMPERATURE: Refrigerate

NEW PREFERRED TRANSPORT TEMPERATURE: Ambient

NOTE: Refrigerate <48 hours is an acceptable transport temperature. Frozen serum is not acceptable.

QUESTIONS: Contact your Mayo Medical Laboratories' Regional Manager or Richard R. Einerson, Mayo Medical Laboratories' Technical Support Telephone: 800-533-1710