



# TEST DEFINITION

9/28/2010

|       |                                 |
|-------|---------------------------------|
| MML   | Name/order code cross-reference |
| Code  | Name                            |
| ----- | -----                           |
| 8413  | Manganese, S                    |

|       |                            |       |  |                   |                         |
|-------|----------------------------|-------|--|-------------------|-------------------------|
| MML   | MML Test setup information |       |  |                   | Perform                 |
| ORDER | EFF                        | TC    | TITLE  | Checking Normals  | Print normals (# coded) |
| CODE  | DATE                       |       |  |                   | Site *                  |
| ----- | -----                      | ----- | -----  | -----             | -----                   |
| 8413  | 8/27/2010                  |       | Manganese, S                                 |                   | SDL                     |
|       |                            |       | Transport temp : Refrig\Ambient OK\Frozen OK |                   |                         |
|       |                            |       | 8413 Manganese, S                            |                   |                         |
|       |                            |       | Units: ng/mL                                 |                   |                         |
|       |                            |       | NO SEX                                       |                   |                         |
|       |                            |       | No age :                                     | ; Not Established |                         |
|       |                            |       | 1D to 15Y :                                  | ; Not Established |                         |
|       |                            |       | 16Y to 150Y : 0.6-2.3                        | ; 0.6-2.3         |                         |
|       |                            |       | MALE   |                   |                         |
|       |                            |       | No age :                                     | ; Not Established |                         |
|       |                            |       | 1D to 15Y :                                  | ; Not Established |                         |
|       |                            |       | 16Y to 150Y : 0.6-2.3                        | ; 0.6-2.3         |                         |
|       |                            |       | FEMALE                                       |                   |                         |
|       |                            |       | No age :                                     | ; Not Established |                         |
|       |                            |       | 1D to 15Y :                                  | ; Not Established |                         |
|       |                            |       | 16Y to 150Y : 0.6-2.3                        | ; 0.6-2.3         |                         |

\*Performing Site Legend

|     |  |               |                                  |
|-----|--|---------------|----------------------------------|
| SDL | Mayo Clinic Dpt of Lab Med & Pathology | LAB DIRECTOR: | Franklin R. Cockerill, III, M.D. |
|     | Superior Dr.                           |               |                                  |
|     | 3050 Superior Dr. NW                   |               |                                  |
|     | Rochester, MN 55901                    |               |                                  |

MML Messages used as normals

|       |       |
|-------|-------|
| CODE  | TEXT  |
| ----- | ----- |

Total of 0 normals codes

\*\*\* End of Report \*\*\*