



TEST DEFINITION

12/6/2011

NAME/ORDER CODE CROSS-REFERENCE
CODE NAME

HTR20 SEROTONIN RECEPTOR 2A/2C, SALIVA

TEST SETUP INFORMATION

ORDER CODE	RESULT CODE	TITLE	CHECKING NORMALS	PRINT NORMALS	PERFORM SITE *
HTR20		SEROTONIN RECEPTOR 2A/2C, SALIVA			MCR
		TRANSPORT TEMP : AMBIENT			
32936		HTR2A -1438G>A	UNITS:		
32937		HTR2A 74C>A	UNITS:		
32938		HTR2A 102 T>C	UNITS:		
32939		HTR2A IVS2 A>G	UNITS:		
32940		HTR2A 1354 C>T	UNITS:		
32941		HTR2C -759 C>T	UNITS:		
32942		HTR2C 796 G>C	UNITS:		
32943		REVIEWED BY	UNITS:		
32944		HTR2 GENOTYPE	UNITS:		

TEST CODE ALWAYS MESSAGE - [NP0118]

NP0118 DIRECT POLYMORPHISM ANALYSIS IS PERFORMED FOLLOWING PCR AMPLIFICATION OF THE TYPE 2A AND 2C SEROTONIN RECEPTOR GENES. ABSENCE OF A DETECTABLE GENE MUTATION OR POLYMORPHISM AT THESE LOCATIONS DOES NOT RULE OUT THE POSSIBILITY THAT A PATIENT HAS A POLYMORPHISM ELSEWHERE IN THE TYPE 2 SEROTONIN RECEPTOR GENES.

TEST CODE ALWAYS MESSAGE - [TLDT]
TLDT LABORATORY DEVELOPED TEST.

*PERFORMING SITE LEGEND

=====
MCR MAYO CLINIC DEPT. OF LAB MED AND PATHOLOGY
 200 FIRST STREET SW
 ROCHESTER, MN 55905

LAB DIRECTOR: FRANKLIN R. COCKERILL, III, M.D.

=====
*** END OF REPORT ***